

# SLE5000 INFANT VENTILATOR

Mode Select **SIMV** Standby

PSV Setup Support  OFF



Apnoea Setup Disabled when backup rate is 20 BPM or higher.

TTV  OFF

Additional pressure support can be given.

No Current Alarms 00

Hist Adjust



II

Print

Light

Ti Meas. 1.00

BPM Tot 21

Trigger 1

Vte (ml) 6.50

Vmin (l) 0.138

Leak (%) 1

Resist. 168.80

Compl. 0.30

C20/C 1.00

Mean P 10



Manual Breath



PROXIMAL AIRWAY

EXHALATION PORT FROM PATIENT

FLOW SENSOR

FRESH GAS TO PATIENT

